

VELLA® Varnish 5% Sodium Fluoride Varnish with NuFluor® (22600ppm Fluoride Ion)

Caution: U.S. Federal Law restricts this device to sale by or on the order of a dentist.

VELLA® Varnish is a 5% sodium fluoride varnish for application to enamel or dentin for the treatment of dental hypersensitivity. 1 mL of **VELLA** contains 50 mg of sodium fluoride (22.8 fluoride ion). **VELLA** Varnish is tolerant of moisture and saliva and will adhere to moist surfaces.

Indications for Use: Dentinal hypersensitivity.

Contraindications: Ulcerative gingivitis, stomatitis.

Precautions: For 24 hours post treatment with **VELLA**, other fluorides such as gels, foams, and rinses should not be administered. The patient should discontinue use of fluoride tablets for 2 – 3 days following treatment.

Warnings: Edematous swelling has been reported only in rare instances, especially following application to extensive surfaces with fluoride varnish. Dyspnea, although extremely rare, has occurred in asthmatic people. Nausea can occur in patients with known sensitive digestive systems, especially following extensive applications.

Administration: For the treatment of dentin hypersensitivity, in most cases, two applications at several days interval are sufficient. The fluoride content of **VELLA** Varnish is such that neither acute nor chronic side effects are to be expected if applied according to instructions.

MSDS: To obtain MSDS visit, www.preventech.com or call Preventech at 1-800-474-8681.

VELLA® & VELLA®Mini Varnish Instructions for Use

1. For optimum results, prior to applying **VELLA** lightly dry excessive saliva by rubbing with a 2 x 2 along the tooth surface. If biofilm or calculus is present, a prophylaxis can be conducted using **nada®** Pumice Paste.
2. For **VELLA**, gently peel open the foil sealing material starting from the unsealed tab on the left side. Lift straight up to avoid tearing. For **VELLAMini**, hold the pouch near the flavor stripe, grab the perforated tear-away strip and slowly pull. Gently squeeze **VELLA** on the back of an exam glove or into a dappen dish.
3. Remove the applicator brush and use the bristle end to mix **VELLA** until it reaches a uniform consistency.
4. Apply a very thin coat of **VELLA** evenly to all of the desired tooth surfaces.
5. Let **VELLA** dry for approximately 20 seconds and then instruct the patient to close their mouth.
6. **VELLA** hardens on contact with saliva so the patient may leave immediately after the application.
7. Instruct the patient to allow **VELLA** to remain in contact with the tooth surface for a minimum of 2 hours. For maximum fluoride uptake, **VELLA** can remain on the treated area overnight and brushed off the next morning.
8. It is recommended that the patient be instructed to only eat soft foods and drink cold liquids and avoid eating hard foods and brushing/flossing during the treatment period.



Manufactured by:
Preventive Technologies, Inc
4330-C Matthews-Indian Trail Road
Indian Trail, NC 28079
1-800-474-8681
704-849-2416
704-849-2417 FAX
www.preventech.com

Made in U.S.A. RM142 Rev C 11/15

VELLA® and **nada®** are registered trademarks of Preventive Technologies, Inc.