

ēsa and the ēsamate prophylaxis system: The reliable, heat sterilizable prophylaxis solution.

ēsa MW-Midwest



ITEM NO.	DESCRIPTION
3300013	Soft Cup—100 per box
3300023	Firm Cup—100 per box
3300033	Tapered Brush—100 per box
3300043	Soft Cup Non-Latex—100 per box
3300047	Soft Cup Non-Latex—1000 per box
3300057	Firm Cup Non-Latex—1000 per box



ēsamate^M_W
Low-Speed Handpiece

ITEM NO. 333334
DESCRIPTION
1 ēsamate^M_W Low Speed Handpiece,
Lube-Free, 5,000 RPM

ēsa ST-Star Titan



ITEM NO.	DESCRIPTION
3300113	Soft Cup—100 per box
3300123	Firm Cup—100 per box
3300133	Tapered Brush—100 per box
3300143	Soft Cup Non-Latex—100 per box
3300117	Soft Cup—1000 per box
3300127	Firm Cup—1000 per box
3300147	Soft Cup Non-Latex—1000 per box



ēsamate^S_T
Low-Speed Handpiece

ITEM NO. 333333
DESCRIPTION
1 ēsamate^S_T Low Speed Handpiece,
Lube-Free, 5,000 RPM

Call your favorite dental dealer to take advantage of the following promotional offers!

**Purchase 5–100 Boxes,
Get 1–100 Box FREE**

**Purchase 1–1000 Box,
Get 1–100 Box FREE**

All ēsa Prophylaxis Cups are now Latex-Free.

**Purchase 1 ēsamate handpiece,
get 1 ēsa MW or ST (100) box FREE**

**Purchase 3 ēsamate^M_W or
ēsamate^S_T Low Speed Handpieces,
Get 1 FREE**

**Purchase 2,000 ēsa MW or ēsa ST,
Get ēsamate HP FREE**


To redeem these “Buy/Get offers” simply place a qualifying order through the dealer of your choice. Then mail, fax or email your dealer invoice to Preventech. Email to customerservice@preventech.com. Allow 4 weeks for delivery. This offer cannot be combined with any other offer. Offer expires 12/31/18.



Your satisfaction is 100% guaranteed. If any Preventech product does not live up to your standards, return it for a full refund. Customer Service: 1-800-474-8681. ēsa is presently only available for Star Titan, Midwest Shorty and Rhino handpieces. All third party marks – ® and ™ – are the property of their respective owners.

PREVENTECH 4330-C Matthews-Indian Trail Road • Indian Trail, NC 28079 • 800.474.8681 • 704.849.2416 • fax: 704.849.2417 • preventech.com

Are you compliant with infection control standards in your state?



- States with specific heat sterilization policy and CE requirements
- States that subscribe to CDC Guidelines and require infection control continuing education
- States that subscribe to CDC Guidelines without requiring infection control continuing education
- States with specific heat sterilization policy for hand piece infection control

ēsa
extended straight attachment

ēsa disposable prophy angles bring super simple compliance to all 50 states.

State Dental Boards are shifting toward sterilization.

CDC Guidelines recommend that, "dental handpieces and the attachments should always be heat sterilized between patients, "adding" if a semicritical item is heat-sensitive, DHCP should replace it with a heat-tolerant or disposable alternative." Additionally, they state, "If a semicritical item does not come in contact with mucous membranes or non-intact skin, it does not require heat sterilization."¹ Most practices utilize a plastic barrier on the handpiece to prevent cross contamination. Currently, 39 state dental boards subscribe to the CDC Guidelines.


In 11 states², state dental practice acts mandate the heat sterilization of handpieces. It's not a recommendation. It's the law. With infection control and the sterilization of instruments becoming more transparent to patients, more states may be on the verge of adopting similar mandates.

With ēsa CDC compliance is only the beginning.

THINK COMPLIANCE – ēsa fits a number of heat sterilization tolerant handpieces: Midwest Shorty, Rhino, Star Titan³ and our ēsamate lube-free handpieces. Using ēsa means you're in compliance no matter where you work.

THINK LIGHTER – ēsa eliminates the nosecone, making it 50% lighter than traditional handpiece connections, and making it easier on your wrist and hands.

THINK EFFICIENCY – Fast and easy changeover between patients. Break free from the drudgery of sterilizing and repairing nose cones.



ēsa MW-Midwest



ēsa ST-Star Titan



ēsamate MW



Midwest® Rhino®



Midwest® Shorty®



ēsamate ST



Star® Titan® 5000 rpm



Star® Titan® 20,000 rpm with torque multiplier

ēsa costs less to use than traditional handpiece and DPA connections.

Per CDC Guidelines	Average DPA	ēsa
Angle Cost	20¢ - 78¢	70¢
Nose Cone Amortization	24¢	0
Plastic Barrier Sleeve	8¢	8¢
Sterilization Pouch	9¢	0
Cleaning Agents	4¢	0
Labor Cost*	85¢	0
Total Cost**	\$1.50 – \$2.08	78¢
Carpal Tunnel Concerns	Possible	Unlikely
100% Guarantee	?	Yes

Source: Henry Schein Dental website (March 2017)

*Cost based on Dental Assistant 2017 average hourly salary of \$16.92. Assumes 3 minutes of infection control preparation. Cost is higher if labor is performed by a Registered Dental Hygienist.

**Does not include handpiece/nose cone repair costs.

Talk to us about the lube-free ēsamate prophy system.

ēsamate handpieces are available in two models: ēsamate ST, with an aluminum housing that weighs just 2 ounces, and our new ēsamate MW, which features a stainless steel housing and weighs 3.2 ounces. Both handpieces are heat sterilization tolerant, offer a 360° swivel, and connect quickly and easily to ēsa disposable prophy angles.

Both handpieces are 5,000 rpm, perfect for prophys, and backed by a 2-year warranty.

THINK COMPLIANCE – The ēsamate handpieces meet CDC Guidelines and all 50 states' dental practice acts! Which means no compliance concerns should your state move from CDC Guidelines to a "heat sterilization only" requirement.

THINK RELIABILITY – Proven air-driven performance. No recharging or batteries!

THINK PEACE OF MIND – Both ēsa DPA and our ēsamate handpieces are backed by our 100% guarantee (if they do not live up to your standards, you may return them for a full refund).

¹Summary of Infection Prevention Practices in Dental Settings, Centers for Disease Control and Prevention, March, 2016, P. 14.

²State Dental Practice Boards in CA, FL, IN, KS, MO, NM, OH, OR, SC, VA and WA require semicritical items to be sterilized after every patient.